

**Supplementary Table S5. Primers used in long-range and multiplex PCR**

<b>Primer</b>	<b>Primer sequence</b>	<b>Reference</b>
Long-range PCR-F	5'-AGTCAGGACACCCCATTTTT-3'	galGal3
Long-range PCR-R	5'-GGTTCTGTATCCACAACAACAGA-3'	galGal3
Primer walk 1F	5'-CGCTGAGATTGAAGCAACAC-3'	
Primer walk 1R	5'-GGAGATCGCGACAGTTTTGT-3'	
Primer walk 2F	5'-TTCCCCTTTTTGAACCACAC-3'	
Primer walk 2R	5'-GGGGTAGCAAGGATGATTGA-3'	
Multiplex PCR F	5'-GCATTCACAAACGGGTGTA-3'	galGal3
Multiplex PCR R1	5'-CAAAACCACAAAGGTAATGTTCA-3'	galGal3
Multiplex PCR R2	5'-CCCAGCAGTAAGCCCTACAT-3'	EAV-HP (NC_005947.1)