

TABLE S3. Probes used in Northern blot analysis

Probe (region recognized)	Sequence (5'-3')
Probe 01 (5'-A ₀)	TCAGGTCTCTGCTGC
Probe 02 (18S)	AGCCATTCGCAGTTTCACTG
Probe 03 (D-A ₂)	TTAAGCGCAGGCCCGGCT
Probe 04 (A ₂ -A ₃)	TGTTACCTCTGGGCC
Probe 05 (5.8S)	GCGTTCTTCATCGATGC
Probe 06 (5'E-C ₂)	TGAGAAGGAAATGACGCT
Probe 07 (E-C ₂)	GGCCAGCAATTTCAAGTTA
Probe 08 (25S)	TACTAAGGCAATCCCGGTTGG
U3	GGATTGCGGACCAAGCTAA
scR1	ATCCCGGCCGCCTCCATCAC