


## RESEARCH ARTICLE

# Sensing chemical-induced genotoxicity and oxidative stress via yeast-based reporter assays using NanoLuc luciferase

Minami Shichinohe<sup>1</sup>, Shun Ohkawa<sup>1</sup>, Yuu Hirose<sup>2</sup>, Toshihiko Eki<sup>1\*</sup>

**1** Molecular Genetics Laboratory, Toyohashi, Japan, **2** Department of Applied Chemistry and Life Science, Laboratory of Genomics and Photobiology, Toyohashi University of Technology, Toyohashi, Aichi, Japan

\* eki@chem.tut.ac.jp



## Abstract

Mutagens and oxidative agents damage biomolecules, such as DNA; therefore, detecting genotoxic and oxidative chemicals is crucial for maintaining human health. To address this, we have developed several types of yeast-based reporter assays designed to detect DNA damage and oxidative stress. This study aimed to develop a novel yeast-based assay using a codon-optimized stable or unstable NanoLuc luciferase (*yNluc* and *yNlucCP*) gene linked to a DNA damage- or oxidative stress-responsive promoter, enabling convenient sensing genotoxicity or oxidative stress, respectively. End-point luciferase assays using yeasts with a chromosomally integrated *RNR3* promoter (<sup>P</sup>*RNR3*)-driven *yNluc* gene exhibited high levels of chemiluminescence via NanoLuc luciferase and higher fold induction by hydroxyurea than a multi-copy plasmid-based assay. Additionally, the integrated reporter system detected genotoxicity caused by four different types of chemicals. Oxidants (hydrogen peroxide, *tert*-butyl hydroperoxide, and menadione) were successfully detected through transient expressions of luciferase activity in real-time luciferase assay using yeasts with a chromosomally integrated *TRX2* promoter (<sup>P</sup>*TRX2*)-linked *yNlucCP* gene. However, the luciferase activity was gradually induced in yeasts with a multi-copy reporter plasmid, and their expression profiles were notably distinct from those observed in chromosomally integrated yeasts. The responses of *yNlucCP* gene against three oxidative chemicals, but not diamide and zinc oxide suspension, were observed using chromosomally integrated reporter yeasts. Given that yeast cells with chromosomally integrated <sup>P</sup>*RNR3*-linked *yNluc* and <sup>P</sup>*TRX2*-linked *yNlucCP* genes express strong chemiluminescence signals and are easily maintained and handled without restrictive nutrient medium, these yeast strains with NanoLuc reporters may prove useful for screening potential genotoxic and oxidative chemicals.

## OPEN ACCESS

**Citation:** Shichinohe M, Ohkawa S, Hirose Y, Eki T (2023) Sensing chemical-induced genotoxicity and oxidative stress via yeast-based reporter assays using NanoLuc luciferase. PLoS ONE 18(11): e0294571. <https://doi.org/10.1371/journal.pone.0294571>

**Editor:** Alvaro Galli, CNR, ITALY

**Received:** August 3, 2023

**Accepted:** November 4, 2023

**Published:** November 22, 2023

**Copyright:** © 2023 Shichinohe et al. This is an open access article distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

**Data Availability Statement:** All relevant data are within the paper and its [Supporting Information](#) files.

**Funding:** The authors received no specific funding for this work.

**Competing interests:** The authors have declared that no competing interests exist.

## Introduction

DNA damage caused by mutagenic chemicals increases the risk of cancer by causing genetic mutations in humans [1]. Oxidative damage is constantly generated by cellular metabolism-

derived reactive oxygen species (ROS) and/or environmental oxidative chemicals [2]. Endogenous ROS and exogenous oxidants damage cellular molecules, such as DNA, leading to cellular dysfunction and genetic mutations; additionally, they are implicated in carcinogenesis, neurodegenerative diseases, cardiovascular diseases, inflammation, and aging [3]. Bioassays have been developed to assess potential genotoxic and oxidative chemicals. To assess potential genotoxicity in eukaryotes, *Saccharomyces cerevisiae*-based reporter assays have been developed using reporter genes linked to several DNA damage-inducible promoters, including *RNR2* [4,5], *RNR3* [5–9], *RAD54* [4,6,10], *RAD51* [11], *HUG1* [9,12], *PLM2*, and *DIN7* [13]. The transcriptional induction of these genes is governed by the activation of the DNA damage checkpoint pathway in response to DNA damage caused by mutagens [14]. The budding yeast *S. cerevisiae* is an eukaryote like humans that can be easily manipulated; therefore, yeast-based reporter assays are suitable for assessing the potential genotoxicity in mammals [15]. However, oxidative damage in animals has been evaluated by measuring the amounts of oxidized biomolecules, such as 8-hydroxy-2-deoxyguanosine, in urine and blood samples [16]. Yu et al. conducted a yeast-based assay for sensing oxidative stress using a redox-sensitive GFP [17]. Several oxidative stress-inducible genes, including *TRX2* and *CTT1* have been identified by genome-wide gene expression analyses in yeast [18–22]. However, only a few studies using yeast-based reporter assays with an oxidative stress-inducible reporter gene have been reported so far [23–26]. Dolz-Edo et al. [23] investigated the transcriptional expression of a destabilized luciferase reporter gene linked to four oxidative stress-responsive promoters (*GRE2*, *CTT1*, *SOD2*, and *CCP1*) in yeast strains following exposure to hydrogen peroxide and menadione. They observed concentration-dependent transient expressions induced by these oxidants. Lewinska et al. [25] and Jayaraman et al. [24] developed yeast reporter strains with *YHB1* and *TRX2* promoter-linked *gfp* genes and observed GFP expression caused by nitric oxide donors and two oxidants (hydrogen peroxide and diamide), respectively.

Previously, we developed a novel yeast-based assay system using both a sensor and an *Escherichia coli lacZ*-reporter plasmid, which could detect low concentrations of genotoxic agents more efficiently than conventional reporter systems [5]. In addition, yeast-based assays for genotoxicity using reporter plasmids carrying the *Cypridina noctiluca* secretory luciferase [8], *GFP*, and firefly luciferase [26] genes linked to the ribonucleotide reductase subunit 3 gene promoter ( $P^{RNR3}$ ) have been successfully used to detect chemical genotoxicity. In addition, we developed a novel yeast-based assay for detecting oxidative agents using a reporter gene encoding unstable firefly luciferase containing the CL1 and proline-glutamate-serine-threonine-rich (PEST) protein destabilizing sequence, linked to the thioredoxin 2 gene promoter ( $P^{TRX2}$ ) [26]. Many oxidative stress-responsive genes are regulated by transcription factors, such as Yap1p [27], and transiently expressed during 100 min after oxidant treatment [23]. Therefore, the system used in our aforementioned study was based on detecting the transient expression of luciferase caused by oxidative agents using the real-time luciferase assay with yeasts and a multi-copy plasmid carrying a  $P^{TRX2}$ -linked destabilized luciferase gene.

However, although genotoxicants and oxidative agents were detected using these plasmid-based reporter assays, one issue remained unresolved. A nutrient-selective medium was required to maintain the plasmids in the yeast during the culture and assay, making the assay time-consuming and inconvenient. In the present study, we generated yeast cells carrying a chromosomally integrated reporter gene using CRISPR/Cas9. The luminescence intensity in the chromosomal integration reporter assay is expected to decrease considerably due to the presence of a single copy of the reporter gene; therefore, a codon-optimized NanoLuc luciferase (yNluc) was used as a reporter to resolve this issue. NanoLuc luciferase is an engineered luciferase originating from the deep-sea shrimp *Oplophorus gracilirostris* [28]. It catalyzes the commercialized substrate Nano-Glo<sup>®</sup> and produces a signal approximately 150-fold brighter

than firefly luciferase in an ATP-independent manner [29]. Masser et al. first reported the development of a yeast-based reporter assay using yNluc and its usefulness [30]; this was followed by the development of a yeast-based estrogen sensor by Cevenini et al. using NanoLuc [31].

Thus, the current study aimed to compensate for the decreased luminescence in the assay using a chromosomally integrated reporter yeast by virtue of its strong luminescence.

## Materials and methods

### Chemicals

Methyl methanesulfonate (MMS), hydroxyurea (HU), mitomycin C (MMC), *tert*-butyl hydroperoxide (*t*-BHP), menadione, diamide, and zinc oxide were purchased from Sigma-Aldrich Inc. (St. Louis, MO). Camptothecin (CPT) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) were obtained from FUJIFILM Wako Pure Chemical Corp. (Osaka, Japan), and phleomycin (Phl) was purchased from InvivoGen (Hong Kong). CPT and menadione were dissolved in dimethyl sulfoxide (Sigma-Aldrich Inc.) and diluted with distilled water. Zinc oxide does not dissolve completely in water; therefore, a sonicated zinc oxide suspension (0.32 g/L) was serially diluted in distilled water.

### Strains

The *S. cerevisiae* strains prepared and used in this study are shown in Table 1. Yeast cells were grown at 30°C in yeast-peptone-dextrose (YPD) medium containing 1% yeast extract, 2% peptone, and 2% glucose [32]. Those with luciferase reporter plasmids were maintained and cultured in synthetic dextrose minimal (SD) medium without histidine (ForMedium, Hunstanton, UK).

### Preparations of yNluc reporter plasmids

Two types of yNluc reporter plasmids (pESC-HISΔGAL1/10<sup>-P</sup>RNR3-yNluc and pESC-HISΔGAL1/10<sup>-P</sup>TRX2-yNlucCP) derived from the multi-copy plasmid vector pESC-HIS (Agilent, Santa Clara, CA) were prepared to detect the genotoxicity and oxidative stress, respectively (Table 2). The yNluc and yNlucCP genes encode a codon-optimized Nluc and an unstable Nluc with the CL1 and PEST protein destabilizing sequence, respectively. Both DNAs were prepared by polymerase chain reaction (PCR) using the KOD FX Neo DNA polymerase (Toyobo, Tokyo), the primer sets (yNluc-5'-F and pESChis2DSma-yNlucend-IFR for yNluc gene, and yNluc-5'-F and pESChis2DSma-yNlucPEST-IFR for yNlucCP gene) (S1 Table), and the yNluc-pUC57-Amp plasmid DNA (Azenta Life Sciences, South Plainfield, NJ) containing a synthesized yNlucCP gene as a template (S1 Fig). The promoter fragments of RNR3 and TRX2 were amplified from yeast genomic DNA (BY4741) using KOD FX Neo, the primer set (pESChis2DSma-RNR3P-IFR and RNR3P-yNluc25-R for RNR3 and pESChis2DSma-TRX2P500-IFR and TRX2P500-yNluc25-R for TRX2), and the promoter-containing plasmid

**Table 1. The yeast strains used in this study.**

Strains	Genotype	Source
BY4741	MATa, his3-Δ1, leu2-Δ0, met15-Δ0, ura3-Δ0	Invitrogen
IMX672	MATa, ura3-52, trp1-289, leu2-3,112, his3Δ, can1Δ::cas9-natNT2	EUROSCRAF
BY4741- <sup>P</sup> RNR3-yNluc	MATa, his3-Δ1, leu2-Δ0, met15-Δ0, ura3-Δ0, can1Δ:: <sup>P</sup> RNR3-yNluc	This study
BY4741- <sup>P</sup> TRX2-yNlucCP	MATa, his3-Δ1, leu2-Δ0, met15-Δ0, ura3-Δ0, can1Δ:: <sup>P</sup> TRX2-yNlucCP	This study

<https://doi.org/10.1371/journal.pone.0294571.t001>

**Table 2. The plasmids used in this study.**

Plasmid	Description
pRS415-LEU-Cas9	A pRS415 plasmid containing <i>cas9</i> expression cassette DNA at <i>Sma</i> I site
pMEL10	An RNA expression plasmid with <i>URA3</i> marker and guide RNA (gRNA) for <i>CAN1</i> gene, which was generated by Mans et al. [33] and obtained from EUROSCARF
<i>yNluc</i> -pUC57-Amp	A pUC57 plasmid containing a codon-optimized <i>yNanolucCP</i> gene, which was synthesized by Azenta Life Sciences (GENEWIZ) based on the nucleotide sequence of pCA955 kindly provided by Prof. Andréasson (Stockholm University) [30]
pESC-HISΔ <i>GAL1/10</i>	A <i>GAL1/10</i> promoter-removed derivative of the pESC-HIS vector [26]
pESC-HISΔ <i>GAL1/10</i> - <sup>P</sup> <i>RNR3-luc2</i>	pESC-HISΔ <i>GAL1/10</i> with a <i>RNR3</i> promoter-linked <i>luc2</i> gene [26]
pESC-HISΔ <i>GAL1/10</i> - <sup>P</sup> <i>TRX2-luc2CP</i>	pESC-HISΔ <i>GAL1/10</i> with a <i>TRX2</i> promoter-linked <i>luc2CP</i> gene [26]
pESC-HISΔ <i>GAL1/10</i> - <sup>P</sup> <i>RNR3-yNluc</i>	A multi-copy <i>yNluc</i> reporter plasmid driven by <i>RNR3</i> promoter, prepared in this study
pESC-HISΔ <i>GAL1/10</i> - <sup>P</sup> <i>TRX2-yNlucCP</i>	A multi-copy unstable <i>yNlucCP</i> reporter plasmid driven by <i>TRX2</i> promoter, prepared in this study

<https://doi.org/10.1371/journal.pone.0294571.t002>

DNA (pESC-HISΔ*GAL1/10*-<sup>P</sup>*RNR3-luc2* and pESC-HISΔ*GAL1/10*-<sup>P</sup>*TRX2-luc2CP* for *RNR3* and *TRX2*, respectively). <sup>P</sup>*RNR3* and <sup>P</sup>*TRX2* fragments were connected to *yNluc* and *yNlucCP* DNAs by PCR using the primer sets (pESChis2DSma-RNR3P-IFF and pESChis2DSma-yNlucend-IFR, and pESChis2DSma-TRX2P500-IFF and pESChis2DSma-yNlucPEST-IFR, respectively) and then cloned into the *Sma* I site of the pESC-HISΔ*GAL1/10* plasmid using the In-Fusion HD cloning kit (Takara, Osaka). The resultant pESC-HISΔ*GAL1/10*-<sup>P</sup>*RNR3-yNluc* and pESC-HISΔ*GAL1/10*-<sup>P</sup>*TRX2-yNlucCP* were used for the genotoxicity and oxidative stress assays, respectively. The nucleotide sequences of the cloned DNAs, along with their flanking regions in the recombinant plasmids, were determined in both directions using a BigDye<sup>®</sup> Terminator v3.1 Cycle Sequencing Kit (Thermo Fisher Scientific Inc., Waltham, MA) and an automated DNA sequencer (Applied Biosystems model 3130xl Genetic Analyzer) by Macro-gen Japan (Tokyo, Japan). The sequence data were assembled and analyzed using the ATGC and Genetyx software (version 13, Genetyx Co., Tokyo). The oligo DNAs used in this study were synthesized by FASMAC (Atsugi, Japan) and Integrated DNA Technologies (Coralville, IO).

### Preparation of chromosomally integrated reporter strains

Two BY4741-derived yeast strains carrying a chromosomally integrated <sup>P</sup>*RNR3-yNluc* and <sup>P</sup>*TRX2-luc2CP* gene (BY4741-<sup>P</sup>*RNR3-yNluc* and BY4741-<sup>P</sup>*TRX2-yNlucCP*) were prepared by replacing *CAN1* with CRISPR/Cas9 (Table 1), as described by Mans et al. [33]. We used CRISPR/Cas9 for chromosomal integration of the reporter construct instead of classical genomic integration because of the high efficiency of integrations at the correct loci and the capability to perform multiple chromosomal integrations without any restrictions regarding the number of nutrient-selective genes. In brief, yeast cells with the Cas9-expression plasmid pRS415-LEU-Cas9 were co-transformed by the *CAN1*-gRNA expression plasmid pMEL10 and the promoter-linked reporter DNA with 5'- and 3'-flanking sequences of *CAN1* at each end. The resultant *LEU*<sup>+</sup>- and *HIS*<sup>+</sup>-transformants were investigated by colony PCR using allele-specific primer sets (S1 Table) with KOD FX Neo polymerase to confirm the integration of the promoter-linked reporter gene at the *CAN1* locus on chromosome V (S2 Table). The DNAs for integration were prepared by amplifying the <sup>P</sup>*RNR3-yNluc* and <sup>P</sup>*TRX2-luc2CP* DNAs from

pESC-HIS $\Delta$ GAL1/10-<sup>P</sup>RNR3-*yNluc* and pESC-HIS $\Delta$ GAL1/10-<sup>P</sup>TRX2-*yNlucCP*, respectively, using KOD FX Neo with the primer set (TADH1-5F2-25CAN1-Rtail and TCYC1-3R3-25CAN1-Ftail). The 5'- and 3'-flanking DNAs were prepared by PCR using KOD FX Neo and the primer sets (5'-CAN1-F and 5-CAN1-R\_5F2\_30tail, and 3-CAN1-F\_3R3\_30tail and 3'-CAN1-R), respectively. Purified promoter-reporter DNA and 5'-flanking DNA were connected by PCR using the primer set (5'-CAN1-F and TCYC1-3R3-25CAN1-Ftail). The resultant DNAs were further connected with 3'-flanking DNA by PCR using the 5'-CAN1-F and 3'-CAN1-R primer set to generate the promoter-linked reporter DNA with the 5'- and 3'-CAN1-flanking sequences.

### Yeast transformation

Yeast cells were transformed with pESC-HIS-derived reporter plasmids using a lithium acetate protocol modified by Gietz et al. [34]. The transformants were selected on SD agar plates without histidine (ForMedium), and independent colonies were streaked onto a fresh selection of agar plates before use. For gene integration, the yeast cells were first transformed by the Cas9-expression plasmid pRS415-LEU-Cas9 with the *LEU2* marker. Then, the *LEU*<sup>+</sup>-transformants were co-transformed by the *CAN1*-gRNA expression plasmid pMEL10 with the *HIS3* marker and a promoter-linked reporter gene fragment as described above. The integrated strains were selected on SD agar plates without leucine and histidine.

### End-point luciferase assay for genotoxicity

Yeast cells with pESC-HIS $\Delta$ GAL1/10-<sup>P</sup>RNR3-*yNluc* were grown on SD agar plates without histidine and cultured with 10 mL of histidine-free SD medium in a 50 mL conical tube at 30°C with continuous shaking for 24 h. Yeast cells with chromosomally integrated <sup>P</sup>RNR3-*yNluc* gene (BY4741-<sup>P</sup>RNR3-*yNluc*) were cultured in 10 mL of the YPD medium. The cells were collected by centrifugation to remove the growth medium and suspended in YPD medium at an absorbance of 600 nm (OD<sub>600</sub>) of approximately 1.0. Subsequently, 100  $\mu$ L of the yeast suspensions were placed in a 96-well white plate (Coster, No. 3912) in triplicates, and the tested chemicals were added at the indicated concentrations. The yeast cells in the microplate were incubated at 30°C for 24 h under saturated humidity in a Tupperware. Next, a 100  $\mu$ L-aliquot of the diluted cell suspension was transferred into a new 96-well white plate containing 10  $\mu$ L of 100-fold-diluted Nano-Glo<sup>®</sup> Luciferase Assay Substrate (Promega, Madison, WI). After 30 min, the A<sub>600</sub> and chemiluminescence intensity were measured using a multimode plate reader (Tecan Infinite M1000, Männedorf, Switzerland). The luciferase activity is shown in arbitrary units defined as the luminescence (counts) in 1 s normalized by the A<sub>600</sub>. The fold induction was calculated as the ratio of the luciferase activity in the presence and absence of each test chemical. In this study, we investigated the responses of a chromosomally integrated <sup>P</sup>RNR3-*yNluc* gene following exposure to four different types of genotoxic chemicals: methyl methanesulfonate (MMS), an alkylating agent, Phleomycin (Phl) that induces DNA double-strand breaks, mitomycin C (MMC) causing DNA strand cross-links, and camptothecin (CPT), a DNA topoisomerase I inhibitor that breaks replication forks.

### Real-time luciferase assay for oxidative stress

Yeast cells with pESC-HIS $\Delta$ GAL1/10-<sup>P</sup>TRX2-*luc2CP* and chromosomally integrated <sup>P</sup>TRX2-*luc2CP* gene (BY4741-<sup>P</sup>TRX2-*luc2CP*) were pre-cultured and cultured in a microtiter plate without chemicals as described in the previous section. A 100  $\mu$ L-aliquot of the diluted cell suspension was transferred into a new 96-well white plate containing 10  $\mu$ L of 100-fold-diluted Nano-Glo<sup>®</sup> Luciferase Assay Substrate, and oxidative chemicals were quickly added at

the indicated concentrations. At time zero, the  $A_{600}$  and chemiluminescence intensity were measured in a multimode plate reader, and the luminescence intensity was subsequently determined every 10 min until 1 h and 20 min at 28°C. The normalized luciferase activity is shown in arbitrary units, as described in the previous section. The relative maximal activity was defined as the percentage of the highest peak luciferase activity at each concentration divided by the highest peak activity among all the concentrations tested. The results were plotted against the concentrations of the oxidative agents.


### Statistical analysis

A two-tailed paired Student's *t*-test using the T.TEST function in Microsoft Excel was employed for the statistical analysis. A *p*-value of less than 0.01 was considered statistically significant.

## Results

### Detection of genotoxicity using yeast strains carrying an *RNR3* promoter-linked *yNluc* gene on a multi-copy plasmid and a chromosome

Two types of reporter yeast strains carrying a multi-copy plasmid with a  $P^{RNR3}$ -linked *yNluc* gene and a chromosomally integrated reporter gene were generated, and the response of the *yNluc* gene in two yeast reporter strains in the presence of hydroxyurea (HU) was investigated. The luminescence intensity in yeasts with the reporter plasmid was more than 10-fold higher than that in the chromosomally integrated reporter yeasts (Fig 1A and 1B, S3 Table). However, significantly higher fold inductions were observed in the chromosomally integrated reporter strain treated with HU than those in the plasmid-based reporter strain (Fig 1C, S4 Table). Based on this observation, we investigated the responses of a chromosomally integrated  $P^{RNR3}$ -*yNluc* gene after exposure to four different types of genotoxic chemicals (MMS, Ph1, MMC, and CPT). The luciferase activities (Fig 2A–2C, S5 Table) and fold inductions (Fig 2E–2G) were increased after treatment with MMS, Ph1, and MMC in a dose-dependent manner.



**Fig 1. Response of  $P^{RNR3}$ -linked *yNluc* gene in two reporter strains after hydroxyurea (HU) treatment.** The luciferase-derived luminescence in yeast cells carrying a multi-copy reporter plasmid (A) and a chromosomally integrated reporter gene (B) were measured after exposure to various concentrations of HU for 24 h. The luminescence intensities were normalized using the  $A_{600}$  value and plotted against the HU concentrations with the standard deviations. The fold inductions in two yeast strains are shown in (C). \* Statistically significant ( $p < 0.01$ , Student's *t*-test). The raw dataset for Fig 1 and statistical analysis are shown in S3 and S4 Tables, respectively.

<https://doi.org/10.1371/journal.pone.0294571.g001>

In particular, yeast cells treated with MMS exhibited 100-fold inductions (Fig 2E). Conversely, CPT treatment did not cause a strong induction in luciferase activity in the reporter yeast cells (Fig 2D).


### Detection of oxidative stress using yeast strains carrying a *TRX2* promoter-linked *yNlucCP* gene on a multi-copy plasmid and a chromosome

We generated two yeast-based reporter assays using <sup>P</sup>*TRX2*-linked unstable *yNluc* gene (*yNlucCP*) on a multi-copy plasmid and a chromosome V (*CAN1* locus) and investigated the response of the reporter gene in two strains treated with three oxidative chemicals: hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), *tert*-butyl hydroperoxide (*t*-BHP), and menadione. The maximal luminescence intensity in yeasts with a reporter plasmid was approximately 4.5-fold higher than the peak activity in a chromosomally integrated reporter yeast in the presence of 0.2 mM H<sub>2</sub>O<sub>2</sub> (Fig 3A and 3B, S6 Table). Interestingly, the resultant induction profiles of the luciferase activity were different in the two reporter strains: the luciferase activity gradually increased in yeasts containing a plasmid with the same reporter gene (Fig 3A), whereas transient peaks of *yNlucCP* were detected in the chromosomally integrated reporter yeasts at 20 min after H<sub>2</sub>O<sub>2</sub> addition (Fig 3B), as observed in our previous study [26]. The highest peak of fold induction was found in yeasts with a chromosomally integrated reporter gene in the presence of 0.2 mM H<sub>2</sub>O<sub>2</sub> (Fig 3C). The “relative maximal activity” of luciferase was defined as the percentage of the highest peak activity at a tested concentration divided by the highest peak activity observed in the experiment in order to determine the induced levels of luciferase against the concentration of the oxidants. As shown in Fig 3D, luciferase induction occurred in a concentration-dependent manner up to 0.2 mM H<sub>2</sub>O<sub>2</sub>; the induction level decreased at a high concentration of 0.4 mM H<sub>2</sub>O<sub>2</sub>.

Similar results were obtained from the assay using *t*-BHP (Fig 4, S7 Table). The luciferase peaks were not observed in the plasmid-based yeast reporter strain up to 100 min after treatment (Fig 4A). However, in the chromosomal integration reporter strain, transient peaks were found at 30–40 min after adding *t*-BHP (Fig 4B). Interestingly, the expression profile of luciferase at the highest concentration of *t*-BHP (3.6 mM) was different from those observed at other concentrations: Despite the presence of a peak around 30–40 min after treatment, the activity demonstrated a continuous increase in the presence of *t*-BHP even after 50 min after the culture period (Fig 4B). The relative maximal activity plateaued at 0.6 mM based on the induction peaks of the expression response caused by *t*-BHP (Fig 4D). Consistent observations were obtained from the assays using menadione-treated reporter yeasts (Fig 5, S8 Table). The expression of luciferase was gradually increased in yeasts with a reporter plasmid (Fig 5A), and transient peaks of luciferase expression were consistently observed in yeasts with an integrated reporter gene (Fig 5B). Luciferase activity measured until 200 min after menadione treatment indicating high levels of activity when the concentration of menadione was high (0.1 mM) in the chromosomally integrated reporter strain (Fig 5B). However, in the plasmid-based reporter assay, the levels of luciferase expression reached a plateau at 70–150 min; the plateau level increased when treated with 0.025 mM menadione and reached its maximum at 0.05 mM, followed by a decrease at 0.1 mM (Fig 5A). Menadione treatment clearly induced the expression of an integrated reporter gene as shown by the fold induction and relative maximal activity in this (Fig 5C and 5D) and another assay (Fig 5E and 5F).

### Response of a chromosomally integrated <sup>P</sup>*TRX2-yNlucCP* gene after treatment with three oxidative chemicals


Finally, we tried to detect the potential for oxidative stress caused by diamide and zinc oxide using yeast cells with an integrated <sup>P</sup>*TRX2-yNlucCP* gene. Although diamide is known to cause



**Fig 2. Nluc activity and fold induction in genotoxicant-treated yeast cells with a chromosomally integrated  $P^{RNR3}\text{-}\gamma\text{Nluc}$  gene.** The luminescence intensity was normalized by the cell density and fold induction in the integrated reporter yeasts after exposure to the indicated concentrations of methyl methanesulfonate (MMS) (A, E), phleomycin (Phl) (B, F), mitomycin C (MMC) (C, G), and camptothecin (CPT) (D, H). The raw dataset for Fig 2 is shown in S5 Table.

<https://doi.org/10.1371/journal.pone.0294571.g002>






**Fig 3. Response of  $P^{TRX2}$ -yNlucCP gene in yeast after exposure to hydrogen peroxide ( $H_2O_2$ ).** The luminescence intensities in yeast cells carrying a multi-copy reporter plasmid (A) and a chromosomally integrated reporter gene (B) were measured in the presence of  $H_2O_2$ ; the normalized intensity is shown with the standard deviations. (C) Fold inductions calculated from the activity of a panel (B) were plotted against the culture period. (D) The relative maximal luciferase activity was calculated as the percentage of each peak activity divided by the highest peak activity and plotted against the concentrations as shown in panel b with the standard deviations. The raw dataset for Fig 3 is shown in [S6 Table](#).

<https://doi.org/10.1371/journal.pone.0294571.g003>

oxidative stress, we failed to detect oxidative stress in our previous study [26]. Zinc oxide nanoparticles have been suggested to cause oxidative stress in animals [35]. In the present study, a typical transient expression of luciferase was found in reporter yeast cells exposed to *t*-BHP as a positive control in luciferase activity (Fig 6A, S9 Table) and fold induction (Fig 6B). The results of the relative maximal activity indicated luciferase induction with *t*-BHP (Fig 6C). However, the luciferase activities in yeast cells exposed to diamide (Fig 6D) and zinc oxide suspension (Fig 6G) were >20-fold lower than that in the *t*-BHP-treated cells. Although weak, luciferase peaks were found even in the absence of test chemicals. No increase in fold induction was observed, and the relative maximal activity did not show any chemical-mediated induction of unstable NanoLuc luciferase with diamide (Fig 6E and 6F) and zinc oxide suspension (Fig 6H and 6I).




**Fig 4. Response of a chromosomally integrated  $P^{TRX2}$ -yNlucCP gene in yeast after exposure to *tert*-butyl hydroperoxide (t-BHP).** The luminescence intensities in yeast cells carrying a multi-copy reporter plasmid (A) and a chromosomally integrated reporter gene (B) were measured in the presence of the indicated concentrations of t-BHP. (C) Fold inductions calculated from the activity of a panel (B) are plotted against the culture period. (D) The relative maximal activities in yeasts treated with the indicated concentrations of t-BHP were calculated from the activity of a panel (B) and plotted against the concentrations with the standard deviations. The raw dataset for Fig 4 is shown in [S7 Table](#).

<https://doi.org/10.1371/journal.pone.0294571.g004>

## Discussion

### Detection of genotoxicity using yeast strains carrying a chromosomally integrated *RNR3* promoter-linked yNluc gene


In this study, we used a codon-optimized NanoLuc luciferase (yNluc) that produces strong luminescence as a reporter in yeast-based assays and a DNA damage-inducible promoter ( $P^{RNR3}$ ) to detect genotoxic chemicals. First, we investigated the response of a reporter gene with HU in two reporter systems:  $P^{RNR3}$ -linked yNluc gene on a multi-copy plasmid and integrated at the *CAN1* locus. As shown in [Fig 1](#), high levels of luminescence intensity were observed in yeasts with a single copy of the reporter gene. Furthermore, a higher fold induction was observed in yeasts with an integrated reporter gene than in yeasts containing a reporter plasmid. While  $10^3$ – $10^5$  counts normalized by  $A_{600}$  were observed as luciferase activity induced by MMS in the assay using yeasts with a chromosomally integrated  $P^{RNR3}$ -firefly luciferase (*luc2*) gene [26],  $10^5$ – $10^8$  order of normalized counts was obtained in the MMS-treated yeast cells carrying an integrated  $P^{RNR3}$ -yNluc gene ([Fig 2](#)). Yeasts with an integrated reporter gene can be cultured and maintained in a non-selective medium; hence, they are suitable for genotoxicity assays compared to plasmid-based reporter assays. The assay using this reporter strain successfully detected four different types of DNA damage caused by chemicals: alkylated base damage (MMS), DNA strand breaks (PhI), DNA strand cross-links (MMC), and replication fork arrest (CPT) ([Fig 2](#)). The fold inductions observed in the assays with MMS



**Fig 5. Response of a chromosomally integrated <sup>P</sup>TRX2-yNlucCP gene in yeast exposed to menadione.** The luminescence intensities in yeast cells carrying a multi-copy reporter plasmid (A) and a chromosomally integrated reporter gene (B) were measured in the presence of the indicated concentrations of menadione. (C) The fold inductions were calculated by the activity of a panel (B) during the culture period. (D) The relative maximal activities in yeasts treated with the indicated concentration of menadione were determined by the activity of a panel (B) and plotted against the concentrations with the standard deviations. (E) and (F) Fold induction and relative maximal activity plots, respectively, obtained from another experiment. The raw dataset for Fig 5 is shown in S8 Table.

<https://doi.org/10.1371/journal.pone.0294571.g005>

and MMC are more than 10- and 2-fold higher than those in our previous study, which used a reporter plasmid with a <sup>P</sup>RNR3-linked secretory luciferase gene [8]; nonetheless, the fold inductions in the system with PhI, HU, and CPT in the current study are almost comparable to the observations of our previous study [8]. Despite testing only five chemicals, our reporter



**Fig 6. Response of a chromosomally integrated  $P^{TRX2}$ -yNlucCP gene in yeast cells treated with three chemicals.** Normalized luciferase activities (A, D, G), fold inductions (B, E, H), and relative maximal activities (C, F, I) in yeasts after exposure to *t*-BHP (A–C), diamide (D–F), and zinc oxide suspension (G–I). The raw dataset for Fig 6 is shown in S9 Table.

<https://doi.org/10.1371/journal.pone.0294571.g006>

assays cover the Ames-positive compounds (Table 3). Regarding the detectable concentrations of these chemicals, the lowest effective concentrations (LECs) in our yNluc reporter assay are comparable or slightly improved compared with the LECs in our previous assays [5,8,26], which used different reporters, as well as the yeast-based GreenScreen assay with a  $P^{RAD54}$ -linked *gfp* reporter gene [36].

These findings suggest that yeast strains carrying an integrated  $P^{RNR3}$ -yNluc gene can be used as an alternative to the currently used yeast-based genotoxicity assays, which employ a

Table 3. Detectable concentrations of genotoxic chemicals compared with our previous systems and other systems.

Promoter-reporter constructs in yeast or a bacterial assay for detecting genotoxicity	Exposure period for endpoint assay	Genotoxic chemicals tested					Reference
		Methylmethane sulfonate (MMS)	Hydroxyurea (HU)	Phleomycin (Phl)	Mitomycin C (MMC)	Camptothecin (CPT)	
<sup>P</sup> RNR3- <i>γNluc</i> gene integrated at the <i>CAN1</i> locus	8 h	0.0025%–0.05% (w/v) <sup>a</sup> (0.0005%–0.05%) <sup>b</sup>	5–50 mM (2.5–50 mM)	0.3–10 μg/mL (0.3–10 μg/mL)	0.01–1 mM (0.01–1 mM); Low FI <sup>c</sup>	5–80 μg/mL (5–100 μg/mL); Low FI	Present study
<sup>P</sup> RNR3- <i>Cypridina luc</i> gene on a single copy plasmid	6, 24 h	0.0025%–0.1% (0.0025%–0.1%) (peak at 0.01% and decrease at >0.01%)	10–100 mM (3–100 mM)	0.3–10 μg/mL (0.3–10 μg/mL)	0.3–1 mM (0.03–1 mM) (at 24 h); Low FI	10–120 μg/mL (10–120 μg/mL) (at 6 h); low FI	Ochi et al. [8]
<sup>P</sup> RNR3- <i>lacZ</i> gene in a multicopy plasmid	6 h	0.0025%–0.1% (0.0025%–0.1%) (peak at 0.02%)	Not tested	Not tested	Not tested	Not detected at 5–160 μg/mL	Ichikawa and Elki [5]
<sup>P</sup> RNR3- <i>luc2</i> gene on a multicopy plasmid	6, 24 h	0.0025%–0.02% (0.001%–0.02%)	Not tested	Not tested	Not tested	Not tested	Suzuki et al. [26]
<sup>P</sup> RAD54- <i>gfp</i> gene in a multicopy plasmid (GreenScreen assay)	16–20 h (overnight)	1.02 μg/mL (0.0001%) (LEC) <sup>d</sup>	593 μg/mL (7.8 mM) (LEC)	12.5 μg/mL (LEC)	200 μg/mL (0.6 mM) (LEC)	Not tested	Cahill et al. [36]
Ames test	–	Yes <sup>e</sup>	Yes	Not found	Yes	Not found	Madia et al. [37]

<sup>a</sup>Range of detectable concentrations

<sup>b</sup>range of tested concentrations

<sup>c</sup>low FI: Low fold induction

<sup>d</sup>LEC: Lowest effective concentrations detected in the GreenScreen assay

<sup>e</sup>Yes: Ames positive.

<https://doi.org/10.1371/journal.pone.0294571.t003>

reporter plasmid. In the future, using DNA repair-deficient mutants as hosts can make this genotoxicity assay more sensitive against low concentrations of genotoxicants, as shown in several studies [9,12,38], including ours [8,26].

### Detection of oxidative chemicals by yeast strains carrying a *TRX2* promoter-linked *γNlucCP* gene

In our previous study, we developed a yeast-based assay using a multi-copy plasmid carrying a <sup>P</sup>*TRX2*-driven unstable firefly luciferase (*luc2CP*) gene and detected oxidative stress via the transient expression of luciferase activity 60–90 min after treatment with the oxidants [26]. However, in the current study, the results of the assay using a multi-copy plasmid with the <sup>P</sup>*TRX2-γNlucCP* gene were different from those reported previously [26]. *Nluc* activities after treatment with oxidants gradually increased during the culture period without transient peaks of activity (Figs 3A, 4A and 5A). The reason for these discrepant observations between assays using a plasmid-based reporter with destabilized firefly and NanoLuc luciferases remains unknown.

Yeast cells carrying a chromosomally integrated <sup>P</sup>*TRX2-γNlucCP* gene effectively sensed three representative oxidants (H<sub>2</sub>O<sub>2</sub>, *t*-BHP, and menadione) but not diamide and zinc oxide (Figs 3–6). These results from four oxidants, except for zinc oxide, are consistent with our previous observations in the assay with a *Luc2CP* reporter [26]. Additionally, the maximal levels of luciferase induction in each assay (5–15 fold induction) are comparable to those in our previous study (3–10) [26]. High levels of luminescence intensity derived from unstable NanoLuc luciferase (*NlucCP*) contribute to the transient expression in assays using yeasts with a single copy of the chromosomally integrated reporter gene. We summarized the detectable concentrations in this reporter system compared with our previous system, which used yeast strains

Table 4. Detectable concentrations of oxidants in three reporter systems.

Oxidant tested	<sup>P</sup> <i>TRX2-yNlucCP</i> gene integrated at the <i>CANI</i> locus	<sup>P</sup> <i>TRX2-luc2CP</i> gene in a multicopy plasmid [26]	<sup>P</sup> <i>CTT1-lucCP</i> gene in a single-copy plasmid [23]
H <sub>2</sub> O <sub>2</sub>	0.05–0.4 mM <sup>a</sup> (0.05–0.4 mM) <sup>b</sup>	0.03–0.1 mM (0.01–0.4 mM)	0.1–0.8 mM (0.025–2 mM)
<i>t</i> -BHP	0.2–3.6 mM (0.2–3.6 mM)	0.03–3.6 mM (0.003–3.6 mM)	Not tested
Menadione	0.025–0.1 mM (0.001–0.1 mM)	0.025–0.1 mM (0.01–1 mM)	0.03–0.17 mM (0.02–0.17 mM)
Diamide	Not detected (0.5–4 mM)	Not detected (0.01–5 mM)	Not tested
Diethyl maleate	Not tested	Not detected (0.1–10 mM)	Not tested
Saturated ZnO suspension	Not detected (1/8–no dilution)	Not tested	Not tested

<sup>a</sup>Range of detectable concentrations

<sup>b</sup>range of tested concentrations.

<https://doi.org/10.1371/journal.pone.0294571.t004>

with a *Luc2CP* reporter on a multicopy plasmid, and the yeast-based assay using a catalase 1 (*CTT1*) promoter-driven unstable luciferase gene developed by other researchers [23] (Table 4). In essence, the detectable oxidants and the range of oxidant concentrations are comparable in our real-time luciferase assays with reporter yeasts. Specifically, hydrogen peroxide, *t*-BHP, and menadione can be detected, but diamide cannot. Additionally, the reporter assay using yeasts with an oxidative stress-responsive *CTT1* promoter-linked luciferase gene was able to detect >100 μM hydrogen peroxide, as well as >30 μM menadione [23]. Furthermore, the recombinant yeast strain BioS-OS1/2, featuring an oxidative stress-inducible GFP reporter, detected 0.3 mM H<sub>2</sub>O<sub>2</sub> [24]. The question arises as to what concentrations of oxidants should be detectable in the assay. Previous studies revealed the 50% effective concentrations (EC<sub>50</sub>) of hydrogen peroxide that caused toxicity (viability, cell death, or apoptosis) in mammalian cells: EC<sub>50</sub> values of 100 and 500 μM in mouse endothelial cell lines [39], <100 μM in rat neuronal B50 cells [40], 1 mM in human neuroblastoma SH-SY5Y cells [41], and 30–500 μM (depending on the exposure duration) in rat C6 glioma cells [42]. It has been reported that adverse effects in neurons were induced by 12–500 μM hydrogen peroxide [43], and 25–50 μM hydrogen peroxide led to apoptosis in human Jurkat T-lymphocytes [44]. Our reporter assays can detect 30–50 μM hydrogen peroxide at the lowest concentrations covering the EC<sub>50</sub> values mentioned above. However, these assays are currently unable to effectively detect lower concentrations of hydrogen peroxide at levels <10 μM. In the future, further improvements in the reporter assays may be necessary, including using yeast strains with disrupted genes involved in cell permeability to enhance sensitivity to oxidants.

In this study, we observed a new finding in the assay for oxidative stress. High levels of luciferase activities increased or tailed after the peak in the assays with *t*-BHP and menadione (Figs 4B and 5B). The oxidant-triggered transcriptional induction of yeast genes, including *TRX2*, which encodes most antioxidants and components of the cellular thiol-reducing pathway, is regulated by a leucine zipper transcription factor Yap1p [27,45,46]. Under oxidative stress-free conditions, Yap1p localizes in the cytoplasm due to rapid transport by the nuclear export receptor Crm1p. Oxidative stress causes redox-dependent conformational changes via intramolecular bond formations of the cysteine residues at the C- and N-terminus regions in Yap1p, thereby avoiding interactions between Yap1p and Crm1p. The resultant nuclear accumulation of Yap1p activates the transcription of the oxidative stress-responsive genes [47]. The transient activation of the *TRX2* promoter by oxidants may be caused by the rapid recovery of the redox state of Yap1p via the overlapped antioxidants and the redox regulation systems. Based on this idea, the phenomena described above were likely caused by a delayed recovery of the redox state of Yap1p due to the excess amounts of oxidants; these phenomena

were commonly observed in the presence of high concentrations of the oxidants (3.6 mM *t*-BHP in Fig 4B, and 0.1 mM menadione in Fig 5B).

Oxidative stress with diamide and zinc oxide was poorly detected in the current study (Fig 6). H<sub>2</sub>O<sub>2</sub>, *t*-BHP, and menadione are producers of ROS *in vivo*, whereas diamide causes oxidative stress via the oxidation of the sulfhydryl groups in glutathione. The different levels of luciferase induction observed in the two groups of oxidants could be their distinct mechanisms of action. Zinc oxide nanoparticles are widely used in industrial applications, including cosmetics; however, recent studies show that zinc oxide nanoparticles induce oxidative stress via ROS production in zebrafish [48] and human hepatic cells [49]. Although oxidative stress was not detected in our yeast-based reporter assay using the zinc oxide suspension, further tests are required to confirm the findings using different types of zinc oxide nanoparticles. The reporter system used in the current study was based on transient expression caused by oxidative stress, which is measured by the real-time reporter assay using yeasts with a stress-responsible promoter-linked destabilized reporter gene. The transient expression of the yeast genes is caused by environmental or chemical stress [23,50]. Therefore, similar yeast-based reporter systems have been developed; for example, for the selective detection of drugs [51]. Using an unstable NanoLuc luciferase as a reporter could improve these yeast-based real-time assays by virtue of the high levels of luminescence intensity.

## Conclusion

In this study, yeast-based assays were developed for detecting chemically-induced genotoxicity and oxidative stress using NanoLuc luciferase as a reporter. Despite the use of a single copy of the reporter gene, the assays using yeasts with chromosomally integrated *RNR3*- and *TRX2* promoter-driven stable and unstable *Nluc* gene successfully assessed the genotoxicity and oxidative stress, respectively, more effectively than the reporter plasmid-based assays. Integrated reporter strains are easily maintained in a non-selective medium; therefore, reporter assays with these strains can be used as a convenient tool for screening potential genotoxic or oxidative chemicals.

## Supporting information

**S1 Fig. Nucleotide sequence of the codon-optimized *yNlucCP* gene.** Nucleotide sequence of the codon-optimized *yNlucCP* gene was shown.

(PDF)

**S1 Table. PCR primers and nucleotide sequences used in this study.**

(PDF)

**S2 Table. Primer sets used for confirmation of chromosomal integration by colony PCR.**

(PDF)

**S3 Table. Raw dataset for Fig 1.**

(PDF)

**S4 Table. Statistical analysis by Student's *t*-test for two reporter systems.**

(PDF)

**S5 Table. Raw dataset for Fig 2.**

(PDF)

**S6 Table. Raw dataset for Fig 3.**

(PDF)

**S7 Table. Raw dataset for Fig 4.**  
(PDF)

**S8 Table. Raw dataset for Fig 5.**  
(PDF)

**S9 Table. Raw dataset for Fig 6.**  
(PDF)

## Acknowledgments

The authors thank Prof. Andréasson (Stockholm University) for providing the nucleotide sequence of pCA955 plasmid, Shogo Matsunobu, Shintaro Tochikawa, Kosuke Yamamoto and other members of the laboratory for technical support and discussions.

## Author Contributions

**Data curation:** Toshihiko Eki.

**Formal analysis:** Minami Shichinohe, Shun Ohkawa, Toshihiko Eki.

**Investigation:** Minami Shichinohe, Shun Ohkawa.

**Methodology:** Minami Shichinohe.

**Project administration:** Yuu Hirose, Toshihiko Eki.

**Supervision:** Toshihiko Eki.

**Visualization:** Toshihiko Eki.

**Writing – original draft:** Toshihiko Eki.

**Writing – review & editing:** Yuu Hirose, Toshihiko Eki.

## References

1. Friedberg EC, Walker GC, Siede W, Wood RD, Schultz RA, Ellenberger T. DNA Repair and Mutagenesis 2nd Edition. Washington, DC.: American Society for Microbiology Press; 2005.
2. Ames BN, Shigenaga MK, Hagen TM. Oxidants, antioxidants, and the degenerative diseases of aging. Proc Natl Acad Sci U S A. 1993; 90(17):7915–22. <https://doi.org/10.1073/pnas.90.17.7915> PMID: 8367443; PubMed Central PMCID: PMC47258.
3. Cooke MS, Evans MD, Dizdaroglu M, Lunec J. Oxidative DNA damage: mechanisms, mutation, and disease. FASEB J. 2003; 17(10):1195–214. <https://doi.org/10.1096/fj.02-0752rev> PMID: 12832285.
4. Afanassiev V, Sefton M, Anantachaiyong T, Barker G, Walmsley R, Wöfl S. Application of yeast cells transformed with GFP expression constructs containing the *RAD54* or *RNR2* promoter as a test for the genotoxic potential of chemical substances. Mutat Res. 2000; 464(2):297–308. [https://doi.org/10.1016/s1383-5718\(99\)00209-0](https://doi.org/10.1016/s1383-5718(99)00209-0) PMID: 10648917.
5. Ichikawa K, Eki T. A novel yeast-based reporter assay system for the sensitive detection of genotoxic agents mediated by a DNA damage-inducible LexA-GAL4 protein. J Biochem. 2006; 139(1):105–12. Epub 2006/01/24. <https://doi.org/10.1093/jb/mvj011> PMID: 16428325.
6. Boronat S, Pina B. Development of *RNR3*- and *RAD54*-GUS reporters for testing genotoxicity in *Saccharomyces cerevisiae*. Anal Bioanal Chem. 2006; 386(6):1625–32. Epub 2006/09/28. <https://doi.org/10.1007/s00216-006-0751-4> PMID: 17004060.
7. Jia X, Zhu Y, Xiao W. A stable and sensitive genotoxic testing system based on DNA damage induced gene expression in *Saccharomyces cerevisiae*. Mutat Res. 2002; 519(1–2):83–92. [https://doi.org/10.1016/s1383-5718\(02\)00129-8](https://doi.org/10.1016/s1383-5718(02)00129-8) PMID: 12160894.
8. Ochi Y, Sugawara H, Iwami M, Tanaka M, Eki T. Sensitive detection of chemical-induced genotoxicity by the *Cypridina* secretory luciferase reporter assay, using DNA repair-deficient strains of *Saccharomyces cerevisiae*. Yeast. 2011; 28(4):265–78. <https://doi.org/10.1002/yea.1837> PMID: 21456053.



9. Wei T, Zhang C, Xu X, Hanna M, Zhang X, Wang Y, et al. Construction and evaluation of two biosensors based on yeast transcriptional response to genotoxic chemicals. *Biosens Bioelectron.* 2013; 44:138–45. <https://doi.org/10.1016/j.bios.2013.01.029> PMID: 23416315.
10. Walmsley RM, Billinton N, Heyer WD. Green fluorescent protein as a reporter for the DNA damage-induced gene *RAD54* in *Saccharomyces cerevisiae*. *Yeast.* 1997; 13(16):1535–45. [https://doi.org/10.1002/\(SICI\)1097-0061\(199712\)13:16<1535::AID-YEA221>3.0.CO;2-2](https://doi.org/10.1002/(SICI)1097-0061(199712)13:16<1535::AID-YEA221>3.0.CO;2-2) PMID: 9509573.
11. Liu X, Kramer JA, Swaffield JC, Hu Y, Chai G, Wilson AG. Development of a highthroughput yeast-based assay for detection of metabolically activated genotoxins. *Mutat Res.* 2008; 653(1–2):63–9. Epub 2008/05/20. <https://doi.org/10.1016/j.mrgentox.2008.03.006> PMID: 18485802.
12. Benton MG, Glasser NR, Palecek SP. Deletion of *MAG1* and *MRE11* enhances the sensitivity of the *Saccharomyces cerevisiae* *HUG1P*-GFP promoter-reporter construct to genotoxicity. *Biosens Bioelectron.* 2008; 24(4):736–41. Epub 2008/08/12. <https://doi.org/10.1016/j.bios.2008.06.033> PMID: 18693109.
13. Bui VN, Nguyen TT, Bettarel Y, Nguyen TH, Pham TL, Hoang TY, et al. Genotoxicity of chemical compounds identification and assessment by yeast cells transformed with GFP reporter constructs regulated by the *PLM2* or *DIN7* promoter. *Int J Toxicol.* 2015; 34(1):31–43. <https://doi.org/10.1177/1091581814566870> PMID: 25691521.
14. Elledge SJ, Zhou Z, Allen JB, Navas TA. DNA damage and cell cycle regulation of ribonucleotide reductase. *Bioessays.* 1993; 15(5):333–9. <https://doi.org/10.1002/bies.950150507> PMID: 8343143.
15. Eki T. Yeast-based genotoxicity tests for assessing DNA alterations and DNA stress responses: a 40-year overview. *Appl Microbiol Biotechnol.* 2018; 102(6):2493–507. Epub 20180208. <https://doi.org/10.1007/s00253-018-8783-1> PMID: 29423630.
16. Kasai H. Analysis of a form of oxidative DNA damage, 8-hydroxy-2'-deoxyguanosine, as a marker of cellular oxidative stress during carcinogenesis. *Mutat Res.* 1997; 387(3):147–63. [https://doi.org/10.1016/s1383-5742\(97\)00035-5](https://doi.org/10.1016/s1383-5742(97)00035-5) PMID: 9439711
17. Yu S, Qin W, Zhuang G, Zhang X, Chen G, Liu W. Monitoring oxidative stress and DNA damage induced by heavy metals in yeast expressing a redox-sensitive green fluorescent protein. *Curr Microbiol.* 2009; 58(5):504–10. <https://doi.org/10.1007/s00284-008-9354-y> PMID: 19184609.
18. Berry DB, Gasch AP. Stress-activated genomic expression changes serve a preparative role for impending stress in yeast. *Mol Biol Cell.* 2008; 19(11):4580–7. <https://doi.org/10.1091/mbc.e07-07-0680> PMID: 18753408; PubMed Central PMCID: PMC2575158.
19. Gasch AP, Werner-Washburne M. The genomics of yeast responses to environmental stress and starvation. *Funct Integr Genomics.* 2002; 2(4–5):181–92. <https://doi.org/10.1007/s10142-002-0058-2> PMID: 12192591.
20. Molina-Navarro MM, Castells-Roca L, Belli G, Garcia-Martinez J, Marin-Navarro J, Moreno J, et al. Comprehensive transcriptional analysis of the oxidative response in yeast. *J Biol Chem.* 2008; 283(26):17908–18. <https://doi.org/10.1074/jbc.M800295200> PMID: 18424442.
21. O'Doherty PJ, Lyons V, Higgins VJ, Rogers PJ, Bailey TD, Wu MJ. Transcriptomic insights into the molecular response of *Saccharomyces cerevisiae* to linoleic acid hydroperoxide. *Free Radic Res.* 2013; 47(12):1054–65. <https://doi.org/10.3109/10715762.2013.849344> PMID: 24074273.
22. Sha W, Martins AM, Laubenbacher R, Mendes P, Shulaev V. The genome-wide early temporal response of *Saccharomyces cerevisiae* to oxidative stress induced by cumene hydroperoxide. *PLoS One.* 2013; 8(9):e74939. <https://doi.org/10.1371/journal.pone.0074939> PMID: 24073228; PubMed Central PMCID: PMC3779239.
23. Dolz-Edo L, Rienzo A, Poveda-Huertes D, Pascual-Ahuir A, Proft M. Deciphering dynamic dose responses of natural promoters and single *cis* elements upon osmotic and oxidative stress in yeast. *Mol Cell Biol.* 2013; 33(11):2228–40. Epub 20130325. <https://doi.org/10.1128/MCB.00240-13> PMID: 23530054; PubMed Central PMCID: PMC3648068.
24. Jayaraman M, Radhika V, Bamne MN, Ramos R, Briggs R, Dhanasekaran DN. Engineered *Saccharomyces cerevisiae* strain BioS-OS1/2, for the detection of oxidative stress. *Biotechnol Prog.* 2005; 21(5):1373–9. <https://doi.org/10.1021/bp050104j> PMID: 16209540.
25. Lewinska A, Grzelak A, Bartosz G. Application of a *YHB1*-GFP reporter to detect nitrosative stress in yeast. *Redox Rep.* 2008; 13(4):161–71. <https://doi.org/10.1179/135100008X259268> PMID: 18647486.
26. Suzuki H, Sakabe T, Hirose Y, Eki T. Development and evaluation of yeast-based GFP and luciferase reporter assays for chemical-induced genotoxicity and oxidative damage. *Appl Microbiol Biotechnol.* 2017; 101(2):659–71. Epub 20161021. <https://doi.org/10.1007/s00253-016-7911-z> PMID: 27766356.
27. Lee J, Godon C, Lagniel G, Spector D, Garin J, Labarre J, et al. Yap1 and Skn7 control two specialized oxidative stress response regulons in yeast. *J Biol Chem.* 1999; 274(23):16040–6. <https://doi.org/10.1074/jbc.274.23.16040> PMID: 10347154.

28. Inouye S, Watanabe K, Nakamura H, Shimomura O. Secretional luciferase of the luminous shrimp *Oplophorus gracilirostris*: cDNA cloning of a novel imidazopyrazinone luciferase. *FEBS Lett.* 2000; 481(1):19–25. [https://doi.org/10.1016/s0014-5793\(00\)01963-3](https://doi.org/10.1016/s0014-5793(00)01963-3) PMID: 10984608.
29. Hall MP, Unch J, Binkowski BF, Valley MP, Butler BL, Wood MG, et al. Engineered luciferase reporter from a deep sea shrimp utilizing a novel imidazopyrazinone substrate. *ACS Chem Biol.* 2012; 7(11):1848–57. Epub 20120830. <https://doi.org/10.1021/cb3002478> PMID: 22894855; PubMed Central PMCID: PMC3501149.
30. Masser AE, Kandasamy G, Kaimal JM, Andréasson C. Luciferase NanoLuc as a reporter for gene expression and protein levels in *Saccharomyces cerevisiae*. *Yeast.* 2016; 33(5):191–200. Epub 20160317. <https://doi.org/10.1002/yea.3155> PMID: 26860732; PubMed Central PMCID: PMC5069653.
31. Cevenini L, Lopreside A, Calabretta MM, D'Elia M, Simoni P, Michelini E, et al. A novel bioluminescent NanoLuc yeast-estrogen screen biosensor (nanoYES) with a compact wireless camera for effect-based detection of endocrine-disrupting chemicals. *Anal Bioanal Chem.* 2018; 410(4):1237–46. Epub 20170930. <https://doi.org/10.1007/s00216-017-0661-7> PMID: 28965124.
32. Dunham M, Gartenberg MR, Brown GW. *Methods in Yeast Genetics and Genomics*. New York: Cold Spring Harbor Laboratory Press; 2015.
33. Mans R, van Rossum HM, Wijsman M, Backx A, Kuijpers NG, van den Broek M, et al. CRISPR/Cas9: a molecular Swiss army knife for simultaneous introduction of multiple genetic modifications in *Saccharomyces cerevisiae*. *FEMS Yeast Res.* 2015; 15(2). Epub 20150304. <https://doi.org/10.1093/femsyr/fov004> PMID: 25743786; PubMed Central PMCID: PMC4399441.
34. Gietz D, St Jean A, Woods RA, Schiestl RH. Improved method for high efficiency transformation of intact yeast cells. *Nucleic Acids Res.* 1992; 20(6):1425. <https://doi.org/10.1093/nar/20.6.1425> PMID: 1561104; PubMed Central PMCID: PMC312198.
35. Sirelkhatim A, Mahmud S, Seeni A, Kaus NHM, Ann LC, Bakhori SKM, et al. Review on Zinc Oxide Nanoparticles: Antibacterial Activity and Toxicity Mechanism. *Nanomicro Lett.* 2015; 7(3):219–42. Epub 20150419. <https://doi.org/10.1007/s40820-015-0040-x> PMID: 30464967; PubMed Central PMCID: PMC6223899.
36. Cahill PA, Knight AW, Billinton N, Barker MG, Walsh L, Keenan PO, et al. The GreenScreen genotoxicity assay: a screening validation programme. *Mutagenesis.* 2004; 19(2):105–19. <https://doi.org/10.1093/mutage/geh015> PMID: 14981157.
37. Madia F, Kirkland D, Morita T, White P, Asturiol D, Corvi R. EURL ECVAM Genotoxicity and Carcinogenicity Database of Substances Eliciting Negative Results in the Ames Test: Construction of the Database. *Mutat Res Genet Toxicol Environ Mutagen.* 2020;854–855:503199. Epub 20200521. <https://doi.org/10.1016/j.mrgentox.2020.503199> PMID: 32660827; PubMed Central PMCID: PMC7374420.
38. Jia X, Xiao W. Compromised DNA repair enhances sensitivity of the yeast *RNR3-lacZ* genotoxicity testing system. *Toxicol Sci.* 2003; 75(1):82–8. Epub 20030612. <https://doi.org/10.1093/toxsci/kgf158> PMID: 12805645.
39. Alamu O, Rado M, Ekpo O, Fisher D. Differential Sensitivity of Two Endothelial Cell Lines to Hydrogen Peroxide Toxicity: Relevance for In Vitro Studies of the Blood-Brain Barrier. *Cells.* 2020; 9(2). Epub 20200210. <https://doi.org/10.3390/cells9020403> PMID: 32050666; PubMed Central PMCID: PMC7072657.
40. Iwata E, Miyazaki I, Asanuma M, Iida A, Ogawa N. Protective effects of nicergoline against hydrogen peroxide toxicity in rat neuronal cell line. *Neurosci Lett.* 1998; 251(1):49–52. [https://doi.org/10.1016/s0304-3940\(98\)00489-3](https://doi.org/10.1016/s0304-3940(98)00489-3) PMID: 9714462.
41. Chetsawang J, Govitrapong P, Chetsawang B. Hydrogen peroxide toxicity induces Ras signaling in human neuroblastoma SH-SY5Y cultured cells. *J Biomed Biotechnol.* 2010; 2010. Epub 20100905. <https://doi.org/10.1155/2010/803815> PMID: 20871828; PubMed Central PMCID: PMC2943129.
42. Gülден M, Jess A, Kammann J, Maser E, Seibert H. Cytotoxic potency of H<sub>2</sub>O<sub>2</sub> in cell cultures: impact of cell concentration and exposure time. *Free Radic Biol Med.* 2010; 49(8):1298–305. Epub 20100729. <https://doi.org/10.1016/j.freeradbiomed.2010.07.015> PMID: 20673847.
43. Hohnholt MC, Blumrich EM, Dringen R. Multiassay analysis of the toxic potential of hydrogen peroxide on cultured neurons. *J Neurosci Res.* 2015; 93(7):1127–37. Epub 20141030. <https://doi.org/10.1002/jnr.23502> PMID: 25354694.
44. Hampton MB, Orrenius S. Dual regulation of caspase activity by hydrogen peroxide: implications for apoptosis. *FEBS Lett.* 1997; 414(3):552–6. [https://doi.org/10.1016/s0014-5793\(97\)01068-5](https://doi.org/10.1016/s0014-5793(97)01068-5) PMID: 9323034.
45. Carmel-Harel O, Stearman R, Gasch AP, Botstein D, Brown PO, Storz G. Role of thioredoxin reductase in the Yap1p-dependent response to oxidative stress in *Saccharomyces cerevisiae*. *Mol Microbiol.* 2001; 39(3):595–605. <https://doi.org/10.1046/j.1365-2958.2001.02255.x> PMID: 11169101.

46. Gasch AP, Spellman PT, Kao CM, Carmel-Harel O, Eisen MB, Storz G, et al. Genomic expression programs in the response of yeast cells to environmental changes. *Mol Biol Cell*. 2000; 11(12):4241–57. <https://doi.org/10.1091/mbc.11.12.4241> PMID: 11102521; PubMed Central PMCID: PMC15070.
47. Kuge S, Arita M, Murayama A, Maeta K, Izawa S, Inoue Y, et al. Regulation of the yeast Yap1p nuclear export signal is mediated by redox signal-induced reversible disulfide bond formation. *Mol Cell Biol*. 2001; 21(18):6139–50. <https://doi.org/10.1128/MCB.21.18.6139-6150.2001> PMID: 11509657; PubMed Central PMCID: PMC87331.
48. Giordo R, Nasrallah GK, Al-Jamal O, Paliogiannis P, Pintus G. Resveratrol Inhibits Oxidative Stress and Prevents Mitochondrial Damage Induced by Zinc Oxide Nanoparticles in Zebrafish (*Danio rerio*). *Int J Mol Sci*. 2020; 21(11). Epub 20200528. <https://doi.org/10.3390/ijms21113838> PMID: 32481628; PubMed Central PMCID: PMC7312482.
49. Pei X, Jiang H, Xu G, Li C, Li D, Tang S. Lethality of Zinc Oxide Nanoparticles Surpasses Conventional Zinc Oxide via Oxidative Stress, Mitochondrial Damage and Calcium Overload: A Comparative Hepatotoxicity Study. *Int J Mol Sci*. 2022; 23(12). Epub 20220616. <https://doi.org/10.3390/ijms23126724> PMID: 35743165; PubMed Central PMCID: PMC9223789.
50. Pascual-Ahuir A, González-Cantó E, Juyoux P, Pable J, Poveda-Huertes D, Saiz-Balbastre S, et al. Dose dependent gene expression is dynamically modulated by the history, physiology and age of yeast cells. *Biochim Biophys Acta Gene Regul Mech*. 2019; 1862(4):457–71. Epub 20190302. <https://doi.org/10.1016/j.bbagr.2019.02.009> PMID: 30836134.
51. Vanadoig-Pedros E, Lozano-Pérez C, Alarcón B, Pascual-Ahuir A, Proft M. Live-cell assays reveal selectivity and sensitivity of the multidrug response in budding yeast. *J Biol Chem*. 2019; 294(35):12933–46. Epub 20190711. <https://doi.org/10.1074/jbc.RA119.009291> PMID: 31296662; PubMed Central PMCID: PMC6721934.