

Table S2

Table S2. Amplification and sequencing primers used.

Amplicon		Primers	
		PCR	Sequencing
A	Fwd	5'-GCTAACCAGCCAGACTCTGGGAGTG-3'	As PCR-primer
	Rev	5'-GGCGAAGGATTTGACTTACAGCAG-3'	As PCR-primer
B	Fwd	5'-GTAGAGATTGAGTGCCTGCTGTTC-3'	As PCR-primer
	Rev	5'-GTGGCTTCCTACCTGAAGCACCG-3'	As PCR-primer
1	Fwd	5'-AAAGAAGACCAGGTTGGAGGCAAAG-3'	5'-GGCAAAGGAGGAAGGGGAGC-3'
	Rev	5'-AGCTAGAGGGCACTTGATGAATAAACG-3'	5'-GGAAATAAGGAAACCACATCAGATGACAC-3'
2	Fwd	5'-TCATATGGACCCTGGCAGTCTCTAATC-3'	5'-TGTGATCCCTTGACTTAACTGTTCTTCC-3'
	Rev	5'-AATCGAGGAAGACAGAATCATCCACC-3'	5'-TCCACCTCTCATTTGCTGCCAG-3'
3	Fwd	5'-AGTGA AAAATGCAAATACTGGTAGGGAGG-3'	5'-CTGGTAGGGAGGGAACACACTAAAATG-3'
	Rev	5'-ACCCAGATGGCCCCTAAAGAGC-3'	5'-CCTCAAAGCTCCAGGACTCTCCAG-3'
4	Fwd	5'-CAAAGGCTGGGATTACACACGAAAC-3'	5'-CAGGATGGAATGAAGACTGGATGC-3'
	Rev	5'-CTGTGTCTGTGCATTTGGATGTGG-3'	5'-CATTTGGATGTGGGGACCAGG-3'
5	Fwd	5'-GGGAACAGGTTTGCCTTCTTGC-3'	5'-TTTGCTTAATGCTGGCCCTTCC-3'
	Rev	5'-GTCTCAGTTTCTAGCATAAGGGGCTG-3'	5'-GGGGCTGTAATCTTCATCTGCAGG-3'
6	Fwd	5'-CCTTCCAGTTAGATCAAACCATTGCTG-3'	5'-CAAACCATTGCTGAAACTGAAGAGGAC-3'
	Rev	5'-CCATCCGCCATTTTGCTGTG-3'	5'-CGCCATTTTGCTGTGTCAATTTCC-3'
7	Fwd	5'-CCTTCTTCCTTTTCCGCCAGG-3'	5'-CCAGGCTTACCATCCAAACAATTCC-3'
	Rev	5'-AGAGCCAGATCACCTTCCACACAC-3'	5'-ACCTCAGGGTGTTGGTTATTCTTTCC-3'
8	Fwd	5'-CTCTTCCAACCTCTGAAAACCATCG-3'	5'-CTCTGAAAACCATCGATGAAGGAATATC-3'
	Rev	5'-TCTGGCTTCCAAACCCTCTTTCTC-3'	5'-CTCTTTCTCTTTGTTTCATCCACCACC-3'
9	Fwd	5'-GTTTTCCCATCTCAGGTGTGTTGC-3'	5'-GGAGACATTGAGGCAGGCACTG-3'
	Rev	5'-GATCTTCCTTGGCCAGGGGTATG-3'	5'-GGGGTATGGGCATGAACAGGAG-3'
10	Fwd	5'-GCATTGTTGGCTGTTCTG-3'	5'-AGTAGAATGGGGGCAGCAC-3'
	Rev	5'-TCCTCACTTGGGACATTTTACAC-3'	5'-GGAGGCTGAAAGCGGAC-3'