

Supplementary Table 2. Coefficient of variation (CV) of ERCC measurements.

Intra-assay; Intra-sample CV	0.20
Intra-assay; Inter-sample CV	0.19

CV of differences is calculated from data presented in Figure 4. $CV = (\text{Standard Deviation of measurements}) / (\text{mean of measurements})$

Intra-assay Intra-sample CV is calculated from the median of intra-assay CV within each sample A-D.

Intra-assay Inter-sample CV is calculated from the median of intra-assay CV across samples A-D.