



Supplemental Figure S1: Comparison of the $\alpha 4\beta 7$ -specific clone Act1 and the $\alpha 4$ -specific clone 7.2R in combination with a separate $\beta 7$ -specific antibody.

A PBMC of healthy controls were stained with either an $\alpha 4$ -specific antibody (clone 7.2R) in combination with a $\beta 7$ -specific antibody (clone FIB504) (left, green) or with the $\alpha 4\beta 7$ -specific antibody, clone Act1 (right, blue). Results were analysed using the Wilcoxon matched-pairs signed rank test, nonparametric test of paired samples. **B** Analysed frequencies of both staining approaches are highly correlated.