



**Figure S5.**

**Cell surface localization of HaPrP molecules expressed in CHO cells.**

CHO cell lines stably expressing wild type and mutant HaPrP were treated with (*lanes 3 and 4*) or without (*lanes 1 and 2*) PI-PLC in order to assess anchorage to the outer leaflet of the plasma membrane via the GPI anchor. Following treatment, proteins in the PI-PLC incubation media (*Media, lanes 1 and 3*) were precipitated, and the cells (*Lysate, lanes 2 and 4*) were harvested in lysis buffer. PrP was detected in the precipitated protein and cell lysate samples by Western blotting.