

Fig. S4: Cytotoxic effect of heat modified citrus pectin at low doses. HepG2 cells were incubated with medium alone (Ctl), 1 μ M staurosporine (STS), 50 μ M etoposide (Etop), different concentrations of hydrolysed citrus pectin (HFCP) or 3 mg/ml citrus pectin (Pectin), with (A) or without (B) 10% feetal calf serum. Cell viability was assessed using a MTT assay after 72h of incubation. Data are means of triplicates +/- SD (n=3). *** : p< 0.001 using ANOVA I and Tukey's multiple comparison test.