

Protein	mutation	efficacy
Akt	overexpression	54.3%
Egfr	activation	100%
	overexpression	100%
Gli	activation	100%
	overexpression	100%
hTert	activation	100%
	overexpression	100%
Nf1	deletion	100%
Nf- κ B	overexpression	2.1%
Pdk1	overexpression	36%
Pi3k	activation	100%
	overexpression	100%
Pkc	activation	100%
	overexpression	100%
Pten	deletion	100%
Ras	activation	100%
	overexpression	100%
Wnt	activation	100%
	overexpression	100%