

Table S4. Residual virus titer (TCID₅₀/mL) after treatment with iota-carrageenan solutions in Diluent P2 (sodium chloride 5 mg/mL adjusted to pH 6-7)

Treatment	Residual virus titer (TCID ₅₀ /mL) after treatment with iota-carrageenan solutions in diluent P2 (sodium chloride 5 mg/mL adjusted to pH 6-7)			
	Replicate 1	Replicate 2	Replicate 3	Median
Untreated	5.62E+05	5.62E+05	5.62E+05	5.62E+05
600 µg/mL	3.16E+01	3.16E+01	3.16E+01	3.16E+01
60 µg/mL	3.16E+01	3.16E+01	3.16E+01	3.16E+01
6 µg/mL	3.16E+01	3.16E+01	3.16E+01	3.16E+01
0.6 µg/mL	3.16E+05	3.16E+04	3.16E+04	3.16E+04
0 µg/mL (Diluent P2)	5.62E+05	3.16E+05	5.62E+05	5.62E+05